

Finger Pulse Oximeter

AOJ-70A

SpO₂ and Pulse Rate can be measured at the same time



Fingertip
Pulse Oximeter



Reminder
Function



Fast
Measurement



LED Screen



SpO₂+PR



Lightweight



Auto Power Off



Accurate
&Reliable



Portable



SpO₂+PR



One machine with two functions

Display Outline



8 Seconds Fast Reading



Fast
Measurement



Different Skin Tones Can Be Measured

Our products have passed clinical trials and widely cover people of different skin tones. The measurement is scientific and accurate.



Lightweight & Portable

Weight: about 43g



36mm

39mm

63mm



SPECIFICATIONS

Mode: AOJ-70A

Display Type: LED

Dimension: 63mm(L)x36mm(W)x39mm(H)

Weight: 43g(Without batteries)

SpO₂ Display Range: 35%-100%

Pulse Rate Measurement Range: 25bpm-250bpm

SpO₂ Accuracy: 70%-100%, $\pm 2\%$;

0%-69%, unspecified

Pulse Rate Accuracy: ± 2 bpm