



Mesh Nebulizer – N200

Features:

- The product will stop working in 10 minutes automatically
- The product has three different nebulization rates for you to choose from.
- The product has a protection mechanism to protect the device from dry-nebulizing.
- The product has a low voltage indication - the indication light will flash when the power is low.

Specifications	
Model	N200
Power Supply	Micro USB or 2*AA batteries
Particle Size	1~5um
Noise	<=50dB
Nebulization Rate	>0.2ml/min (Normal/Soft/Strong)
Medication Capacity	8ml

Packing List:

- 1 x device
- 1 x adult mask
- 1 x child mask
- 1 x mouthpiece (optional)
- 1 x usb cable
- 1 x user manual