

Mesh Nebulizer – N400



Features:

- The product will stop working in 10 minutes automatically
- The product has a built-in rechargeable battery set.
- The product has a protection mechanism to protect the device from dry-nebulizing.
- The product has a low voltage indication - the indication light will flash when the power is low.

Specifications

Model	N400
Power Supply	Micro USB or 2*AA batteries
Particle Size	1~5um
Noise	<=50dB
Nebulization Rate	>0.2ml/min
Medication Capacity	8ml

Packing List:

- 1 x device
- 1 x adult mask
- 1 x child mask
- 1 x mouthpiece (optional)
- 1 x usb cable
- 1 x user manual