

Electronic Blood Pressure Monitor

Model: AOJ-30F



2x99
Memories




Irregular
Heartbeat
Detection



Audio & Visual
Warning



Keep Still
Indication



Accurate
Measurement



Detection for
Cuff Wearing




Display Date
and Time



Large Screen
Display

Large four colors backlit LCD display




White for Standby



Green for Normal



Orange for Mild High Pressure or Hypotension



Red for High Pressure

Technical parameters

Measurement Method: Oscillometric

Measurement Range: Pressure 20-295 mmHg (0-39.33kPa)
Pulse rate 40-199 heart beats/min

Measuring Accuracy: Pressure: ± 3 mm Hg (± 0.4 kPa) ;
Pulse: $\pm 5\%$ of reading

Power supply: 4xAAA, DC6V

Low power parameters: $5.2V \pm 0.1V$ prompt low power,
 $< 4.8V \pm 0.1V$ shutdown

Blood pressure test current:
 < 30 mA (when the air pump valve is at non-working state)

Shutdown current: < 50 uA

Automatic shutdown time: 1 minute

Display : Large four colors backlit LCD display

Cuff size: 22-40cm

Unit size: 118mmX98mmX62.7mm

LCD size: 65mmX50mm (3.2 inches)

Weight: 300g (battery not included)