

Compressor Nebulizer – YS-03



Features:

- Adopted high-efficient oil-free piston pump
- Small and round size, portable
- Ideal supplies for family and medical institution
- Low noise, environment protection items
- Using molecular filters to fully filter the air pollutants

Specifications

Model	YS-03
Power Supply	220V/50Hz
Particle Size	1~5um
Noise	<65dB
Nebulization Rate	>0.2ml/min
Medication Capacity	2-12ml

Packing List:

- 1 x device
- 1 x mask
- 1 x medication cup
- 1 x air tube
- 1 x user manual