

2022 Best Seller AOJ-35D


# Wrist Blood Pressure Monitor

Higher Accuracy with 9 Excellent Performances.


To Check Your Blood Pressure and Heartbeat Anytime Anywhere.



# EXCELLENT PERFORMANCES



Irregular  
Heartbeat  
Detection



Good-performance  
Chips



Voice  
Broadcast



Super Large  
LCD Screen



One-touch  
Operation



Portable and  
Lightweight



Two Users Mode  
(2\*99 Memories)



4-color  
Backlit Display



Auto-power  
Off in 20s

# One Button Operation

Easy-to-use with one-button operation.  
The whole process just needs less than  
1 minute.



# Really Portable Really Lightweight

It weighs only 128g with sleek and stylish wearable design, making it really portable and really comfortable to carry when going out.

128g

128g (without batteries)



about half size of your  
mobile phone



# Large LCD Screen with 4-color Backlit Display


With big fonts, the screen will change color according to the measured blood pressure value. Easy to read even in darkness and friendly for all users.



- White for standby
- Green for normal
- Yellow for mild high blood pressure of hypertension
- Red for high blood pressure

# Technology-driven Quality with Higher Accuracy

Good-performance chips and upgraded algorithm make measurement results more accurate.



# Voice Broadcast

Voice broadcast during the whole measurement process.

Suitable for the seniors and people with poor vision. The voice function can be turned off if you don't need it.

Your blood pressure is  
Systolic pressure 109 mmHg,  
Diastolic pressure 78 mmHg.

Your pulse is 88 bpm.



# Dual User Mode

## 2X99 Sets of Memories

The BP monitor can store 99 readings for each user and each reading comes with date and time stamp for easy checking.



One device  
for 2 persons



Dual user  
switching





# Intelligent Inflation Technology

Intelligent inflation technology is used to avoid discomfort caused by the sudden change of pressure and ensures a comfortable measurement.



# Auto-power Off Function

It will shut down automatically without use for 20 seconds. Save power and prolong the device's lifespan.



**Energy-efficient:  
Auto-powers off in  
20 seconds**


# Adjustable Cuff with Skin-friendly Materials

**Long adjustable cuff with skin-friendly materials**

Fit adult arms from 5.31 IN to 7.68 IN (13.5 cm~19.5 cm) in circumference



**Environmentally friendly ABS plastic shell**



# PARTS IDENTIFICATION



- ① User Selection/General Settings Button
- ② Memory Button
- ③ HD LCD Screen
- ④ START/STOP Button
- ⑤ Wrist Cuff

