

Blood Glucose Monitor, Accurate Testing.

G-425-2



Fast

Exact

Stable

Four electrode technique, Anti-interference, More Accurate

Better interference resistance of RBC accretion avoid the influence for testing result, make the measurement result more accurate.



Less blood sample, fast and stable testing

0.1ul blood sample.
only need 5 second from testing to evaluation.



Automatic stripping to avoid cross infection

After the test, press the left button,

One button automatic stripper, no manual contact, convenient and hygienic




Considerate design for better using

448 Set memory

No Coding

Auto strip Ejection

Small size, easy to carry



448 sets Memory for easily observing Blood glucose trending

No Coding
Insert the strip, meter will recognize the strip automatically to avoid misdetecting caused by coding.



Auto strip Ejection
Press one button to eject strip automatically, no need to collect used strip. Convenient and sanitary

Small and portable, easy to carry and make test anytime and anywhere.



Basic Parameter

Storage Environment	-20°C ~55°C, RH ≤93%
Memory	448sets
Testing unit	mmol/L or mg/dL
Data Reading	LCD screen
Testing time	5S
Product size	90 x 54 x 13cm
Weight	73g (including battery)
Testing range	1.1 ~33.3mmol/L
Targe blood	Microvascular blood
Battery	2 x AAA Battery, DC3V

How to Use

