

# BLOOD GLUCOSE TEST STRIPS

HIGH ACCURACY | NO CODING | AUTO EJECTION



The first Glucose Test Strip manufacturer in Guangdong Province.



1 High Accuracy

2 Automatically absorb blood

3 0.7ul blood sampling

4 No coding

## Product information



### Basic Parameter

Storage Environment : 4°C ~30°C, RH≤80%

Testing time: 5S

Packing: 25pcs/vial ; 2 vials / box

Testing range: 1.1 ~33.3mmol/L

01

Biosensing technique, only 5S testing time.



Four electrode technique, Anti-interference, More Accurate  
Better interference resistance of RBC accretion avoid the influence for testing result, make the measurement result more accurate.

02



03

Automatically absorb blood sample, easy and convenient.



04

批号: B200102F26  
生产日期: 2020.xx  
有效日期: 2022.xx

Each box of strip has batch number, exp date, product date. 2 years warranty.

